


STL File Export Instructions

EXPORTING THE STL FILES FROM ITERO OR INTRA ORAL SCANNER

- Scan on iRecord setting if possible
- Capture the vestibule and palate
- Make sure that you have captured the bite scans in correct occlusion. An accurate wax up and planning cannot be done if the bite is not captured correctly. Review scan to check occlusion.
- Make sure that there are no holes or distortion on the teeth or Palate. Holes or distortion may result in a guide that does not fit. Review the scan for holes/missing information.
- AI should be turned off to capture posterior ridge for missing 2nd molars. Capture the entire ridge of the edentulous area. Manually trim cheek tissue from obscuring the crowns



HOLES ON MODEL



CAPTURE ENTIRE RIDGE OF EDENTULOUS AREA AS FAR BACK AS POSSIBLE

ITERO SCANS - SHARE CASE DIRECTLY THROUGH ITERO MACHINE

- iTero Lab ID (C&M Surgical Guides) : 414424 manufacturing@cdental.com
- Make a note on the form submission that you shared the IOS through Itero directly
- Make sure to attach zipped DICOM to [STL upload form](#) see upload instructions on page 3

STL File Export Instructions

EXPORTING THE STL FILES FROM INTRA ORAL SCANNER

3SHAPE SCANS: SHARE CASE DIRECTLY THROUGH TRIOS MACHINE

- 3shape share Trios files to C&M Surgical Guides, inc manufacturing@cdental.com

OR, TO EXPORT STLS YOURSELF:

If the iTero was scanned as restorative, choose Implants and Other exports-Open Model

Download Settings 1 Scan(s) x

Restorative (1)

CAD/CAM System

Implants and Other Exports - Open Model ▼

Hide patient name

Cancel Download

Choose these options if available

Export Type Open Shell ▼

Data Format File per Arch (arches oriented in occl... ▼

File Type STL+PLY ▼

Hide patient name

Cancel Export

STL File Export Instructions


EXPORTING THE STL FILES FROM INTRA ORAL SCANNER

UPLOADING THE DATA FILES TO CDENTAL WEBSITE

By uploading the case you are agreeing to service fees for case preparation of your scans for surgical planning. Fee lists are available upon request.

STEP 1

GO TO [CASE UPLOAD PAGE ON CDENTAL WEBSITE](#)



STEP 2

WAX UPS ARE ENCOURAGED FOR ACCURATE PLANNING. WAX UPS ARE \$25 PER TOOTH. BE SURE TO INCLUDE WAX UP SITES FOR MISSING TEETH IF NEEDED

STEP 3


CHOOSE ONE OF THE FOLLOWING OPTIONS



STEP 4


BROWSE TO THE DOWNLOAD FOLDER AND ATTACH THE ORTHOCAD EXPORT FOLDER THAT YOU EXPORTED FROM WWW.MYITERO.COM

See page 1 for export instructions.



STEP 5

MAKE SURE TO SELECT THE ONE WITH THE DATE AND TIME THAT YOU EXPORTED THE FILE



STEP 6

Make sure to HIT SUBMIT



If you have any questions, please contact our Sara Tarazi - Chief Technologist

✉ starazi@cdental.com

☎ 510-368-7240