

CDental/MDI Export Instructions – DICOM Only


EXPORTING THE DICOM FILES FROM A CBCT SCAN TO BE MERGED WITH ITERO TAKEN BY C-DENTAL/MDI



Please contact your CBCT machine sales representative to learn how to export DICOM files in the following format. Here is an example for Carestream.


STEP 1

CHOOSE SETTINGS AS SLICE THICKNESS .300 OR .400, SINGLE ARCH



STEP 2

EXPORT WILL BE MULTIPLE .DCM FILES (HUNDREDS OF FILES)



STEP 3


THE .DCM FILES SHOULD BE CONTAINED INSIDE A FOLDER WITH THE PATIENT'S NAME

STEP 4

PLEASE DO NOT INCLUDE MULTIPLE FOLDERS OR OTHER NON DICOM CASE FILES

STEP 5

ZIP THE PATIENT FOLDER RIGHT CLICK, SEND TO COMPRESSED (ZIPPED) FOLDER



CDental/MDI Export Instructions – DICOM Only

EXPORTING THE DICOM FILES FROM A CBCT SCAN TO BE MERGED WITH ITERO TAKEN BY C-DENTAL/MDI CONT.

UPLOADING THE DICOM TO C-DENTAL/MDI FOR MERGE WITH A INTRA ORAL SCAN TAKEN AT C-DENTAL OR MDI

STEP 1

FILL OUT THE UPLOAD FORM ON THE WWW.CDENTAL.COM WEBSITE

[Click Here for Upload Page](#)

STEP 2

WAX UPS ARE ENCOURAGED FOR ACCURATE PLANNING. WAX UPS ARE \$25 PER TOOTH. BE SURE TO INCLUDE WAX UP SITES FOR MISSING TEETH IF NEEDED

STEP 3


CHOOSE ONE OF THE FOLLOWING OPTIONS

WAX-UP : Virtual Wax-Up - CBCT + Intraoral Scan (indicated to...

NON WAX-UP : Merge DICOM & STL Only

STEP 4

BROWSE TO WHERE YOU SAVED THE DICOM FOLDER AND ATTACH THE ZIPPED PATIENT FOLDER WITH THE DICOM

 PATIENT NAME DICOM FOLDER

Attachment - Please zip DICOM files into one folder and use the patient name as the title of the folder**

Drag and drop files here

Browse

STEP 5

Make sure to HIT SUBMIT




Please note on the RX when referring the patient to us for intra oral scan that you will be sending the DICOM. Patients will be responsible for the planning fee which is incorporated into the intra oral scan fee at the time of imaging. ***If we are not aware that you will be sending the DICOM we will need to bill you for the planning fee.***

If no imaging was done by C-Dental/MDI, please submit your case through the C&M website at www.cmsurgicalguides.com/send-us-a-case. C&M cases are for outside imaging and there are instructions for STL + DICOM upload.



If you have any questions, please contact our Sara Tarazi - Chief Technologist

 starazi@cdental.com

 510-368-7240